


A.M.I.[®]- dock TissueBag DLU* 10 mm

- For safe removal of specimen during laparoscopic procedures.
- Reduces waste.



A.M.I.[®]- dock FasciaCutter DLU 10 mm

For safe enlargement of trocar incisions, to accommodate removal of specimen.



*DLU: Disposable Loading Unit

The **A.M.I.[®]- dock System** continues to expand:

The new Loading Units!


C::J MEDICAL.....

Agency for Medical Innovations

A.M.I.®- dock

The Loading Units:


Handy Loading-Units - sterile packed - 20 pieces per Box




Blunt dissector DLU



Pleuroscratcher DLU




Easy-Loop DLU




Endo-Curette DLU




Spatula DLU



Endo-Brush DLU



Scale DLU




TissueBag DLU

NEW!



FasciaCutter DLU

NEW!



Order Code:

5 mm A05 251 Ideal for blunt dissection during minimal invasive procedures:
10 mm A10 251 Securely fastened and stays in shape during usage.

5 mm A05 321 Easy and controlled traumatization of parietal pleura

5 mm A05 451 Tissue bundles and vessels can be ligated under enhanced control due to loop fixation.

10 mm A10 221 To remove and obtain tissue samples, especially plaque tissue.

5 mm A05 721 For precise, controlled mobilization of thin tissue layers.
10 mm A10 181

5 mm A05 361 To obtain cytological probes

5 mm A05 301 For measurement and documentation purposes.

10 mm A10 211 For safe removal of specimen during laparoscopic procedures.
10 units / box

10 mm A10 191 For safe enlargement of trocar incisions, to accommodate removal of specimen.

Examples of indications:

Lap. Cholecystectomy
Lap. Nissen Fundoplication
Thoracoscopic procedures
Diagnostic procedures

Pleurodesis

Lap. Cholecystectomy
Lap. Appendectomy
Colon Resection
Adhesiolyses
Lap. Lymphnode-Dissection

Lap. Acute Abdomen
Thoracoscopic procedures
Pleurectomy
Diagnostic Laparoscopy

Lap. Inguinal Hernia Repair
Lap. Nissen Fundoplication
Lap. Adhesiolyses
Thoracoscopic procedures

Diagnostic Laparoscopy

Lap. Inguinal Hernia Repair
Lap. Nissen Fundoplication
Lap. Varicocele

Lap. Cholecystectomy
Lap. Appendectomy
Lap. Myomectomy
Lap. Oophorectomy

Lap. Cholecystectomy
Lap. Appendectomy

A.M.I.®- dock Applicator Short 5 mm / 10 mm

Applicator short; for pediatric surgery

A.M.I.®- dock Applicator Thoracic 5 mm / 10 mm

Applicator curved; for thoracic surgery

A.M.I.®- dock Applicator 5 mm / 10 mm

The universal instrument in minimal invasive surgery

CJ MEDICAL

CJ Medical

12 Long Wall, Haddenham
Buckinghamshire HP17 8DL, UK
Phone and Fax: +44 (1844) 299501
E-Mail: info@cjmedical.com
www.cjmedical.com