


A.M.I.® Gastric Band Cutter



**Facilitates gastroscopic
removal of migrated Gastric
Band Implants!**

- no surgical intervention
- minimally invasive
- safe, fast and cost-effective

C J MEDICAL.....

Agency for Medical Innovations

A.M.I.® Gastric Band Cutter

A.M.I.® Gastric Band Cutter

Gastric Band Implants are put around the stomach for the treatment of Morbid Obesity.


3 - 5 % of these implanted Gastric Bands do migrate into the stomach and have to be removed.

The A.M.I.® Gastric Band Cutter System facilitates the removal of migrated Gastric Band Implants from the stomach by means of gastroscopy:


- easy,
- fast,
- safe, and
- cost-effective.



A.M.I.® Gastric Band Cutter System
(Gastric Band not included!)



Picture showing Handle and Toggle put in place for cutting.



Stainless Steel Cutting Wire for single use.

The system consist of:

- GBC-Handle
- GBC-Toggle
- GBC-Toggle Holder
- GBC-Flexible Wire Guide
- GBC-Cutting Wires
- Using Manual

Order Code: GBC 360

Stainless Steel Cutting Wire for single use
with the GBC 360
Box of 3 Units

Order Code: GBC 364

CJ MEDICAL

CJ Medical

12 Long Wall, Haddenham
Buckinghamshire HP17 8DL, UK
Phone and Fax: +44 (1844) 299501
E-Mail: info@cjmedical.com
www.cjmedical.com