




PUSEN USCOPE 6.3 Fr SINGLE-USE WITH SUCTION DISS™

DIRECT-IN-SCOPE-SUCTION





- ① **DIRECT-IN-SCOPE-SUCTION** PROVIDES SURGEON CONTROL
- ② **BUILT IN TWO-WAY TAP** INCLUDED WITH EVERY STERILE DEVICE
- ③ **STRAIGHT WORKING CHANNEL** FACILITATES EASY INSTRUMENT ACCESS
- ④ **285° DUAL ANGULATION** ENABLES POWERFUL DEFLECTION
- ⑤ **DUAL LIGHT SOURCE** ELIMINATES SHADOWING AND BLIND SPOTS
- ⑥ **HYDROPHILIC COATED SHAFT**
 - **OUTSIDE DIAMETER 6.3 Fr** WORKING CHANNEL 3.6 Fr
 - **1:1 TORQUE** PROVIDES SUPERIOR SURGEON MANOEUVRABILITY
 - **HD IMAGE PROCESSOR** PV300 PLUG & PLAY RECORD STILL IMAGES AND VIDEO FOOTAGE AND CONNECT TO HOSPITAL STACK SYSTEMS
 - **AVAILABLE IN REVERSE DEFLECTION & STANDARD DEFLECTION**



Product Code	Product Description	Diameter (Fr)	Working Channel (Fr)	Box Quantity (Units)
PU630A-S	Uscope single-use 6.3fr DISS Standard deflection	6.3	3.6	10
PU630A-R	Uscope single-use 6.3fr DISS Reverse deflection	6.3	3.6	10

CUSTOMER SERVICES

ORDER DIRECT: OFFICE@CJMEDICAL.COM

CJ Medical

Charles House, Three Milestone, Truro, Cornwall, TR4 9FB

020 3198 2420

 [@cjmedicaluk](https://twitter.com/cjmedicaluk)

 [cjmedical.com](https://www.cjmedical.com)

 [CJ Medical UK](#)