





# PUSEN USCOPE 6.3 Fr SINGLE-USE WITH SUCTION DISS™

DIRECT-IN-SCOPE-SUCTION





- 1 **DIRECT-IN-SCOPE-SUCTION** PROVIDES SURGEON CONTROL
- 2 **BUILT IN TWO-WAY TAP** INCLUDED WITH EVERY STERILE DEVICE
- 3 **STRAIGHT WORKING CHANNEL** FACILITATES EASY INSTRUMENT ACCESS
- 4 **285° DUAL ANGULATION** ENABLES POWERFUL DEFLECTION
- 5 **DUAL LIGHT SOURCE** ELIMINATES SHADOWING AND BLIND SPOTS
- 6 **HYDROPHILIC COATED SHAFT**
- **OUTSIDE DIAMETER 6.3 Fr** WORKING CHANNEL 3.6 Fr
- **1:1 TORQUE** PROVIDES SUPERIOR SURGEON MANOEUVRABILITY
- **HD IMAGE PROCESSOR** PV300 PLUG & PLAY RECORD STILL IMAGES AND VIDEO FOOTAGE AND CONNECT TO HOSPITAL STACK SYSTEMS
- **AVAILABLE IN REVERSE DEFLECTION & STANDARD DEFLECTION**



Product Code	Product Description	Diameter (Fr)	Working Channel (Fr)	Box Quantity (Units)
PU630A-S	Uscope single-use 6.3fr DISS Standard deflection	6.3	3.6	10
PU630A-R	Uscope single-use 7.5fr DISS Reverse deflection	6.3	3.6	10

## CUSTOMER SERVICES

ORDER DIRECT: [OFFICE@CJMEDICAL.COM](mailto:OFFICE@CJMEDICAL.COM)

### CJ Medical

Charles House, Three Milestone, Truro, Cornwall, TR4 9FB

020 3198 2420

 @cjmedicaluk

 [cjmedical.com](http://cjmedical.com)

 CJ Medical UK