




# PUSEN USCOPE 7.5 Fr SINGLE-USE WITH SUCTION DISS™

DIRECT-IN-SCOPE-SUCTION





- 1 **DIRECT-IN-SCOPE-SUCTION** PROVIDES SURGEON CONTROL
- 2 **BUILT IN TWO-WAY TAP** INCLUDED WITH EVERY STERILE DEVICE
- 3 **STRAIGHT WORKING CHANNEL** FACILITATES EASY INSTRUMENT ACCESS
- 4 **285° DUAL ANGULATION**
- 5 **DUAL LIGHT SOURCE** ELIMINATES SHADOWING AND BLIND SPOTS
- 6 **HYDROPHILIC COATED SHAFT**
  - **OUTSIDE DIAMETER 7.5 Fr** WORKING CHANNEL 3.6 Fr
  - **1:1 TORQUE** PROVIDES SUPERIOR SURGEON MANOEUVRABILITY
  - **HD IMAGE PROCESSOR** PV300 PLUG & PLAY RECORD STILL IMAGES AND VIDEO FOOTAGE AND CONNECT TO HOSPITAL STACK SYSTEMS
  - **PROMOTES GETTING IT RIGHT FIRST TIME (GIRFT)**
  - **AVAILABLE IN REVERSE DEFLECTION & STANDARD DEFLECTION**



## Customer Services

Order Code: PU3033AH-REVERSE  
 Description: Pusen 7.5 Fr DISS Ureteroscope RD with Suction  
 Box Quantity: 5 units

PU3033AH-STANDARD  
 Pusen 7.5 Fr DISS Ureteroscope SD with suction  
 5 Units

Order Direct: office@cjmedical.com  
 Order via Framework: NHS Supply Chain

**CJ Medical**  
 Charles House, Three Milestone, Truro, Cornwall, TR4 9FB  
 020 3198 2420

@cjmedicaluk  
 cjmedical.com  
 CJ Medical UK