

HALO LIFT





HALO™ LIFT

-  **HALO™- HALO™ Identification Enables Precise Visual Guidance For Accurate Patient Positioning**
-  **Single-Valve System Provides Automated Sequential Inflation**
-  **Patient Straps For Patient Security & Stability During Positioning & Lifting**
-  **8-Point Positioning Handles Promote Improved Control**
-  **450Kg Max Weight**
-  **PTFE Low-Friction Base**
-  **Variable Speed Air Blower With Intuitive Dial For Precise Control Engineered For Extended Reliability With 5 Year Service Schedule**

Product Codes	Description	Box Quantity (Units)
HALOLIFT32	HALO™ Air Patient Lift 32"	1
HALOLIFT39	HALO™ Air Patient Lift 39"	1
HALOVSAB	HALO™ Variable Speed Air Blower	1

Customer Services

Order Direct: office@cjmedical.com

CJ Medical
Charles House, Three Milestone, Truro, Cornwall, TR4 9FB
020 3198 2420

 [cjmedicaluk](#)

 [cjmedical.com](#)

 [CJ Medical UK](#)