




PUSEN CYSTOSCOPE SINGLE-USE WITH SUCTION

*See clearly, Diagnose Accurately,
Operate Confidently*





- ① **BUILT IN TWO-WAY TAP** INCLUDED WITH EVERY STERILE DEVICE
- ② **STRAIGHT WORKING CHANNEL** FACILITATES EASY INSTRUMENT ACCESS
- ③ **SUCTION CONTROL** FOR IMPROVED VISIBILITY & SURGEON CONTROL
- ④ **210° UP / DOWN DEFLECTION** FACILITATES **RETROGRADE VIEWS**
- ⑤ **DUAL LIGHT SOURCE** ELIMINATES SHADOWING AND BLIND SPOTS
 - **MINI PCNL APPROVED** COMPATIBLE WITH MINI PCNL SHEATH
 - **EFFICIENT STENT REMOVAL**
 - **1:1 TORQUE** PROVIDES SUPERIOR SURGEON MANOEUVRABILITY
 - **REDUCED OUTER DIAMETER** PROVIDES IMPROVED PATIENT COMFORT
 - **HD IMAGE PROCESSOR** PV300 PLUG & PLAY RECORD STILL IMAGES & VIDEO FOOTAGE CONNECT DIRECT TO HOSPITAL STACK SYSTEM
 - **PROMOTES GETTING IT RIGHT FIRST TIME (GIRFT)**



Customer Services

Order Code: PC200-AR
 Description: Cystoscope Single-Use
 Box Quantity: 10 units

Order Direct: office@cjmedical.com
 Order via Framework: NHS Supply Chain

CJ Medical
 Charles House, Three Milestone, Truro, Cornwall, TR4 9FB
 020 3198 2420

 @cjmedicaluk

 cjmedical.com

 CJ Medical UK