

# activeCut



## The prostate morcellator

Experience the new AMX morcellator for anatomic endoscopic enucleation of the prostate. Designed with surgeons and procedures in mind, it combines excellent technology with an outstanding user experience.



treatment-  
induced UI



one-click collection  
system



integrated  
aspiration



integrated  
scale

The prostate morcellator activeCut enables minimally invasive treatment of **benign prostate obstruction** like never before. Thanks to its intuitive graphical user interface, simple aspiration set-up and integrated scale, the morcellator offers **unique usability** in the operating theatre.

The device features uninterrupted treatment without the need to empty a vacuum reservoir and a **reduced risk of blockage** due to fewer plug connections. While its atraumatic blade geometry ensures patientsafety. The **Active Medical Xcellence GmbH** is an German company, developing practical, secure and intuitive solutions for urology.

## Technical specifications

| activeCut (console)                      |  | REF  | 100-0000                     | 1 |
|------------------------------------------|--|------|------------------------------|---|
| device model                             |  |      | table top device             |   |
| aspiration                               |  |      | integrated aspiration pump   |   |
| classification (acc. MDR)                |  |      | IIa                          |   |
| display                                  |  |      | colored 10" LCD with touch   |   |
| electrical requirements                  |  |      | 100-240 Vac, 50/60 Hz, 90 VA |   |
| activeDrive (handpiece)                  |  | REF  | 100-0002                     | 2 |
| rotation speed                           |  |      | 500–3000 rpm                 |   |
| operation mode                           |  |      | rotation/oscillation         |   |
| sterilisation method                     |  |      | moist heat (acc. ISO 17665)  |   |
| cable length                             |  |      | 3 m                          |   |
| MaCu Operation Kit 10 (consumables)      |  | REF  | 400-0003                     |   |
| <b>activeBlade</b> (blade)               |  | 1 x  |                              | 3 |
| dimensions                               |  |      | O 15 Fr x 385 mm free length |   |
| uses                                     |  |      | reusable (up to 10 x)        |   |
| sterilisation method                     |  |      | moist heat (acc. ISO 17665)  |   |
| <b>activeCollect</b> (collection system) |  | 10 x |                              | 4 |
| capacity                                 |  |      | 300 g                        |   |
| uses                                     |  |      | single use, non sterile      |   |



### CUSTOMER SERVICES

ORDER DIRECT: [OFFICE@CJMEDICAL.COM](mailto:OFFICE@CJMEDICAL.COM)

**CJ Medical**

Charles House, Three Milestone, Truro, Cornwall, TR4 9FB  
**020 3198 2420**

 @cjmedicaluk

 [cjmedical.com](http://cjmedical.com)

 CJ Medical UK