


Equipped with ViewHub
S13 Video Processor

Compatible with both GorgeView rhinolaryngoscope and LoopView bronchoscope, adapting for different scenarios, including outpatient clinics, bedside, ICU, and operating rooms.



Mode	R22	R30	R48
Insertion Portion			
Insertion tube O.D.	2.2 mm	3.0 mm	4.8 mm
Channel width	/	/	2.2 mm
Working length	300 mm	300 mm	365 mm
Angulation range	145 ° U/D	145 ° U/D	210 ° U/D
Optical System			
Resolution	160 K	920 K	920 K
Field of View	120 °	140 °	140 °
Depth of Field	3-100 mm	3-100 mm	3-100 mm
Direction of View	0 °	0 °	0 °

CUSTOMER SERVICES

ORDER DIRECT: OFFICE@CJMEDICAL.COM

CJ Medical
Charles House, Three Milestone, Truro, Cornwall, TR4 9FB
020 3198 2420

@cjmedicaluk
cjmedical.com
CJ Medical UK

MACROLUX



GorgeView
Single-Use Digital Flexible
Rhinolaryngoscope
R22/R30/R48

○ Gentle Access




Pistol handle design for comfortable experience



2.2mm Ultra Slim



Big angle
small radius



More natural drop
proves better softness*

○ Anytime, Anywhere


 Support DICOM **/USB
/HDMI signal

 No cross infection


 Always accessible




○ Vivid Image




Cyan



Chromvue



Limpid




Amber

 Smart Shot**
One click, ten photos
AI choose the best

 One click for image
180° reverse

* Internal insertion tubedrop comparison with traditional scope.

** This feature is expected to be updated with the processor software by the end of 2025.



Diagnostic

R22
Ultra slim

R30
HD

Operation

R48 with
2.0mm
working channel

○ Total Solution