




PUSEN USCOPE 5.1 Fr WORKING CHANNEL DISSTM



PU400



- ❶ **5.1 FR WORKING CHANNEL FACILITATES IMPROVED VISABILITY & ASPIRATION**
- ❷ **150% GREATER INFLOW**
- ❸ **BULLET TIP FOR IMPROVED INSERTION**
- ❹ **INTEGRATED T CONNECTOR**
- ❺ **DISS™ DIRECT IN SCOPE SUCTION**
- ❻ **270° UP / DOWN DEFLECTION**
- ❼ **STRAIGHT WORKING CHANNEL FACILITATES EASY INSTRUMENT ACCESS**
- ❽ **DUAL LIGHT SOURCE ELIMINATES SHADOWING AND BLIND SPOTS**
- ❾ **HD IMAGE PROCESSOR PLUG & PLAY SYSTEM, COMPATIBLE WITH HOSPITAL STACK SYSTEM**
- ❿ **1:1 TORQUE PROVIDES SUPERIOR SURGEON MANOEUVRABILITY**



Product Code	Product Description	Working Length (mm)	Box Quantity (Units)
PU400A	Single-Use 9.2 Fr DISS™ Standard Deflection	650	10
PU400A	Single-Use 9.2 Fr DISS™ Reverse Deflection	680	10

CUSTOMER SERVICES

ORDER DIRECT: OFFICE@CJMEDICAL.COM

CJ Medical

Charles House, Three Milestone, Truro, Cornwall, TR4 9FB

020 3198 2420

 [cjmedicaluk](https://twitter.com/cjmedicaluk)

 [cjmedical.com](https://www.cjmedical.com)

 [CJ Medical UK](https://www.linkedin.com/company/cj-medical-uk)