

# **HALO™** PATIENT FLAT LIFT KIT





1

## **HALO™ LIFT**



**HALO™– HALO™ Identification Enables Precise Visual Guidance For Accurate Patient Positioning**



**Single-Valve System Provides Automated Sequential Inflation**



**Patient Straps For Patient Security & Stability During Positioning & Lifting**



**8-Point Positioning Handles Promote Improved Control**



**450Kg Max Weight**



**PTFE Low-Friction Base**

2

## **HALO™ MAT**



**HALO™– HALO™ Identification Enables Precise Visual Guidance For Accurate Patient Positioning**



**Patient Straps For Patient Security & Stability During Positioning & Lateral Transfer**



**8-Point Positioning Handles Promote Improved Control**



**450Kg Max Weight**

### 3 HALO™ VSAB



Variable Speed Air Blower With Intuitive Dial For Precise Control  
Engineered For Extended Reliability With 5 Year Service Schedule

### 4 HALO™ PPP



The Portable Power Pack Provides Cordless Mobility For Use  
Anywhere Without Mains Power

### 5 HALO™ FLK BAG & TROLLEY



The HALO™ Flat Lift Kit Bag Ensures Smooth, Effortless  
Maneuverability For Efficient Transportation

### 6 HALO™ CART



The HALO™ Cart Provides Safe & Stable Transport, Allowing  
Caregivers To Move Equipment Quickly & Efficiently

| Product Codes | Description                   | Box Quantity (Units) |
|---------------|-------------------------------|----------------------|
| HALOFLK32     | HALO Patient Flat Lift Kit 32 | 1                    |
| HALOFLK39     | HALO Patient Flat Lift Kit 39 | 1                    |

## Customer Services

Order Direct: [office@cjmedical.com](mailto:office@cjmedical.com)

**CJ Medical**  
Charles House, Three Milestone, Truro, Cornwall, TR4 9FB  
**020 3198 2420**



[cjmedicaluk](#)



[cjmedical.com](http://cjmedical.com)



[CJ Medical UK](#)