



Innovators . Educators . Facilitators

Patient Flat Lift Kit Empowered Choice

The Most Efficient and Safest Method for Recovering Fallen Patients from Floor Level with a Suspected Spinal or Pelvic Injury or Post-Resuscitation Events. The Evacuation HoverJack Patient FLK is also available as a speciality Stairwell and Rescue Device



[Tap to watch Animation Video](#)

Patient Flat Lift Kit

Features	Benefits
No weight restrictions	Safe for all Patients to include plus-sized Patients
Lift In Supine Position	Meets with UK Resuscitation guidelines preventing further patient injury during any lift
CPR	CPR can be performed on the FLK
Two Operators	Minimum of Two Caregivers Required to Operate the FLK
Safe Operation	Dignified and safe Patient Lifting and Transfer
Promotes Dignity	Minimised Risk of injury to Caregivers and Patients
Rescue	Can be used within a Rescue Setting
Extra Equipment Not Required	Patients can be rescued from a Bathroom or confined spaces where a hoist will not access or operate
Interchangeable Modular Design	Storage Boxes can be interchanged & configured to Suit
Portable Power Pack (PPP)	Battery Powered
Improved Ergonomic Design	PPP can be added to the FLK Cart at a future date

NHS Supply Chain Codes

Product Name	Product Description	NHS Supply Chain Code
Patient Flat Lift Kit 32 Powered	Transfer Equipment Patient Flat Lift Kit 32" Portable Battery Powered	TLC1312
Patient Flat Lift Kit 39 Powered	Transfer Equipment Patient Flat Lift Kit 39" Portable Battery Powered	TLC1306
Flat Lift Kit 32 Evacuation Powered	Transfer Equipment Patient Flat Lift Evacuation Kit Portable Battery Powered	TLC1309
Flat Lift Kit 32	Transfer Equipment Patient Flat Lift Kit 32"	TLC321
Flat Lift Kit 39	Transfer Equipment Patient Flat Lift Kit 39"	TLC581
Patient Flat Lift Kit Evacuation	Transfer Equipment Patient Flat Lift Evacuation Kit	TLC939
Powered Storage Cart	Transfer Equipment Storage Cart & Portable Battery Powered	TLC1311

Available via Frameworks or Direct from CJ Medical for Volume Discounts



Innovators . Educators . Facilitators

Customer Services

T: 020 3198 2420 E: office@cjmedical.com

CJ Medical, Charles House, Three Mile Stone, Truro, Cornwall, TR4 9FB